

# TRACHEO S.E.T. Safe & Easy Tracheostomy

Una pinza innovativa nella tracheostomia



Una pinza che ha accoppiato la metodica Griggs e quella di Ciaglia, traendo da entrambi i pregi eliminandone i difetti



## CARATTERISTICHE

- Apertura pinza per compressione e non trazione
- Presenza di una valva che permette il passaggio del filo guida senza mai perderlo
- Conformazione pinza analoga a un monodilatatore



## VANTAGGI

- Miglior controllo nel realizzare lo stoma
- Possibilità di regolare la massima apertura delle pinze
- Procedura rapida e semplice per l'inserimento della cannula

# TRACHEO S.E.T. Safe & Easy Tracheostomy

The innovative forceps in Tracheostomy



The dilating forceps that allows the combination of the Griggs and the Ciaglia technique, gaining their full advantages and eliminating their disadvantages



## FEATURES

- The forceps works through compression and not traction
- A valve allows the passage of the guide wire, without losing it
- The shape of the forceps is similar to a mono-dilator



## ADVANTAGES

- Better control during the stoma formation
- Possibility to adjust the maximum opening of the forceps
- Fast and simple procedure for cannula insertion