

USB Video-Rhino-Laryngoscope

RSX / RSX-P Paediatrics



Video-Rhino-Laryngoscope

Model:	RSX	RSX-P
Diameters:	3.9mm	2.9mm
Working length:	310mm	
FOV:	90°	
Tip Articulation:	130°/130°	
Weight:	approx. 330g incl. USB-Connection cable	
Interface:	USB	
Lightsources:	Integrated high power LED lightsource	
	Completely immersible machine cleaning/ disinfection	



Medical Tablet

Model:	10,1" Medical Tablet PC, fanless design
Capacity:	256GB SSD
Operating system:	Windows 10
Power supply:	Exchangeable battery
Resolution:	1280x800 Pixel
Operation:	Multi-Touch-Panel reinforced glass front
Interfaces:	W-LAN 802.11 a/b/g/n, Bluetooth, USB
Speakers:	integrated
Microphone:	integrated
Dimensions:	292x196x20mm
Weight:	1.1Kg


USB Video-Rhino-Laryngoscope



RSX / RSX-P Paediatrics

Plug & Play: All you need

Efficient Diagnostics. Made in Germany




- RSX 3.9mm diameter
- RSX-P 2.9mm diameter (Paediatrics)
- ultra-compact with integrated LED
- USB connection




Configurations




Portable system



directly connection to Medical PC



directly connection to Medical Tablet



Medical Tablet